

MEDI-EVAC MAT ^{CE}

(Part of the “**SKI**” range of products)

The Medi-Evac Mat has been designed to be used as an evacuation aid for a non-ambulant person in an emergency. Similar in design to that of the popular Bariatric Evacmat, the Medi-Evac Mat’s tough, slippery Polyethylene base board allows it to be moved easily over indoor and outdoor surfaces for both horizontal and vertical evacuation.

Features:

- Securely cocoons the user for a greater sense of security while being moved – especially useful with users who may display challenging behaviour while being evacuated due to dementia or mental illness.
- Tough base board enables the Medi-Evac Mat to be moved over any surface – indoor or outdoor.
- Adjustable length to accommodate users of varying heights.
- Three patient securing straps with heavy duty plastic buckles.
- Head and foot end pulling handles.
- Multiple positioning handles
- Four anchor points.
- Internal pillow
- Self-contained in its own storage bag.
- Wall mountable.
- Easy clean after use.



Specifications

Packed: 51cm x 25cm x 24cm approx

Unrolled: 145cm x 240cm approx

Base: Polyethylene board

Material: Flame retardant anti-bacterial PVC

Weight: 5kg

Working load limit: 160kg



HOSPITAL AIDS
Evacuation Equipment Specialists

1 Kingsthorpe Business Centre, Tel: +44(0)1604 586501
63 Studland Road, Kingsthorpe, Fax: +44(0)1604 586576
Northampton, NN2 6NE. England sales@hospitalaids.co.uk
Registered in England No 1777125 www.hospitalaids.co.uk



MADE IN THE UK
DISTRIBUTED WORLDWIDE